


# Impedance : Other Transfer Functions (TF)




## Background on EIS


### Other Transfer Functions:

- Photo-electrochemical TF (IMVS/IMPS)
- Electro-hydrodynamic TF (EHD)
- Magneto-hydrodynamic TF
- Electro-gravimetric TF
- Raman spectro TF
- Thermo-electrochemical TF (TEC)
- Pneumato-chemical Impedance Spectroscopy (PIS)
- ...

### Summary



Current response has the same frequency with an **amplitude  $\delta I$**  and **phase  $\Phi$**   
 Perturbation in potential, (it is also possible to perform the same in Galvano)



Increasing the frequency  moving away from the steady state  $I_{ss}$  vs E curve

$$Z(f) = \mathcal{L}E(t)/\mathcal{L}I(t)$$


$\mathcal{L}$ : Laplace Transform

The impedance is a complex number:

$$Z = a + jb = \text{Re}(Z) + j\text{Im}(Z) \text{ (with } j^2 = -1)$$


$$Z = \rho(\cos\varphi + j\sin\varphi) \text{ with } \rho \text{ the modulus and } \varphi \text{ the phase}$$

## Nyquist diagram



In the Nyquist plot, the impedance for each frequency is plotted in the complex plane  $-\text{Im}(Z)$  vs  $\text{Re}(Z)$ .

## Bode diagram




In the Bode Plot, the modulus and the phase of the impedance are plotted against the frequency of the modulation.

## What if we replace U and/or I by other input/output quantities? Not only electrical quantities.

It is not only electrical Impedance Spectroscopy but Transmittance spectroscopy  
This allows one to study the dependence (output) of the system to an input

It is a measurement of magnitude and phase of the output as a function of frequency, in comparison to the input.

Unit is not in Ohm but is the *ratio* of the two quantities  
 $H = \text{output}/\text{input}$




## What if we replace U and/or I by other input/output quantities? Not only electrical quantities.

It is not only electrical Impedance Spectroscopy but Transmittance spectroscopy  
This allows one to study the dependence (output) of the system to an input

It is a measurement of magnitude and phase of the output as a function of frequency, in comparison to the input.


Unit is not in Ohm but is the *ratio* of the two quantities  
 $H = \text{output}/\text{input}$



$$F(\text{input}) = \text{input} \times \text{TF} = \text{input} \times \frac{\text{output}}{\text{input}} = \text{output}$$

What about if we replace E and/or I by other input/output quantities?

Not only electrical quantities.



Potential  
Current  
Temperature  
Pressure  
Rotation speed of the electrode  
Magnetic field  
Light intensity  
...

Potential  
Current of the electrode  
Current of the secondary electrode  
Temperature  
Pressure  
Magnetic field  
...

## Other Transfer Functions:

	Input	Output	ref
<b>IMPS/IMVS</b>	Light intensity/ $\phi$	PhotoCurrent/I Photovoltage/E	Halme, J. Phys. Chem. Chem. Phys., 2011, 13, 12435–12446
<b>Electro-hydrodynamic TF (EHD)</b>	Rotation speed/ $\Omega$	Current/I Voltage/E	Tribollet <i>et al.</i> J. of Electroanal. Chem. , 2004, Vol. 572, 2, Pages 389–398
<b>Magneto-hydrodynamic TF</b>	Magnetic field/B	Current/I Voltage/E	Olivier <i>at al.</i> J. of Electrochem. Soc., 2004, Vol 151, 2, C112-C118,
<b>Electro-gravimetric TF</b>	Quartz resonator frequency/F or adsorbed mass/m Sauerbrey relationship: $\Delta F = k\Delta m$	Voltage/E	Perrot <i>et al.</i> The Journal of Physical Chemistry B, 2002, Vol. 106, 12, 3182-3191
<b>Raman</b>	Raman intensity/Count.s <sup>-1</sup>	Current/I Voltage/E	Deslouis <i>et al.</i> Electrochimica acta, 2010, 55, 6299-6307
<b>Thermo-electrochemical TF (TEC)</b>	Temperature/T	Current/I Voltage/E	Olivier, A. <i>et al.</i> Electrochimica Acta, 1996, Vol 41, 17, 2731-2736
<b>Pneumato-chemical Impedance Spectroscopy (PIS)</b>	Pressure/P	Voltage/E	Millet, P. <i>et al.</i> J. Phys Chem B; 2005 109 24016-24024




## Other Transfer Functions:

	Input	Output	ref
IMPS/IMVS	Light intensity/ $\phi$	PhotoCurrent/I Photovoltage/E	Halme, J. Phys. Chem. Chem. Phys., 2011, 13, 12435–12446
Electro-hydrodynamic TF (EHD)	Rotation speed/ $\Omega$	Current/I Voltage/E	Tribollet <i>et al.</i> J. of Electroanal. Chem. , 2004, Vol. 572, 2, Pages 389–398
Magneto-hydrodynamic TF	Magnetic field/B	Current/I Voltage/E	Olivier <i>at al.</i> J. of Electrochem. Soc., 2004, Vol 151, 2, C112-C118,
Electro-gravimetric TF	Quartz resonator frequency/F or adsorbed mass/m Sauerbrey relationship: $\Delta F = k\Delta m$	Voltage/E	Perrot <i>et al.</i> The Journal of Physical Chemistry B, 2002, Vol. 106, 12, 3182-3191
Raman	Raman intensity/Count.s <sup>-1</sup>	Current/I Voltage/E	Deslouis <i>et al.</i> Electrochimica acta, 2010, 55, 6299-6307
Thermo-electrochemical TF (TEC)	Temperature/T	Current/I Voltage/E	Olivier, A. <i>et al.</i> Electrochimica Acta, 1996, Vol 41, 17, 2731-2736
Pneumato-chemical Impedance Spectroscopy (PIS)	Pressure/P	Voltage/E	Millet, P. <i>et al.</i> J. Phys Chem B; 2005 109 24016-24024

Each type of Transfer Functions is a disciplin itself.  
We will see only the two first

# IMVS/IMPS



$\phi$ : light intensity

$$\text{IMVS} = \mathcal{L}\Phi(t) / \mathcal{L}E(t)$$



$$\text{IMPS} = \mathcal{L}\Phi(t) / \mathcal{L}I(t)$$

I: PhotoCurrent - IMPS

U: PhotoVoltage - IMVS

**IMVS:** Intensity **M**odulated photo**V**oltage **S**pectroscopy

**IMPS:** Intensity **M**odulated **P**hoto**C**urrent **S**pectroscopy




## Which system?

Solar cell, Grätzel cell,


## Which information?

- $D \approx a \cdot d^2 / \tau_{IMPS}$   
where  $a \approx 0.393$  for weakly absorbed modulated light  
D: diffusion coefficient of electrons/ $m^2 \cdot s^{-1}$   
d: photoelectrode film thickness/m
- $\tau_{IMPS} \approx \tau$
- $L = (D \tau)^{1/2}$   
L: electron diffusion length/m
- $\eta_{COL} \approx 1 - (\tau_{IMVS} / \tau_{IMPS})$   
 $\eta_{COL}$ : Efficiency electron collection



$$IMVS = \mathcal{L}\Phi(t) / \mathcal{L}E(t)$$

$$IMPS = \mathcal{L}\Phi(t) / \mathcal{L}I(t)$$




$$IMVS = \mathcal{L}\Phi(t) / \mathcal{L}E(t)$$



$$IMPS = \mathcal{L}\Phi(t) / \mathcal{L}I(t)$$

In that case, there is a relationship between EIS, IMPS and IMVS.

For the same condition (same E, same I, same light intensity),

$$EIS = IMVS / IMPS$$






Lifetime of the electron  $\tau_n$  is related to the characteristic frequency  $f_c$  by the following equation

$$\tau_n = 1/(2 \pi f_c)$$


$$f_{\min} = 1/(2 \pi RC)$$

Resistance, Capacitor, Minimum frequency and Electron lifetim of white LED vs. light intensity.

Light intensity/ W.m <sup>-2</sup>	R/ Ohm	C/ μF	f <sub>min</sub> / Hz	τ <sub>n</sub> / ms
0.500	67.55	84.04	28.05	5.94
0.750	52.45	90.48	33.56	4.59
1.000	42.12	99.91	35.32	4.04
1.250	32.84	110.78	43.77	3.55
1.500	24.73	124.73	51.62	3.12
1.750	18.95	146.48	57.36	2.74
2.000	14.27	161.35	69.20	2.41
2.250	11.09	172.00	83.52	1.87
2.500	8.83	185.32	97.29	1.64
2.749	7.22	193.07	114.30	1.45
2.999	6.14	201.40	128.69	1.27



# Electro-HydroDynamic TF (EHD)



$\Omega$ : rotation speed

$$EHD_E = \mathcal{L}E(t) / \mathcal{L}\Omega(t)$$


$$EHD_I = \mathcal{L}I(t) / \mathcal{L}\Omega(t)$$

I: current

U: voltage

$$Z = EHD_E / EHD_I$$





## Which system?

Useful for analysis electrochemical system that are either partially or completely limited by mass transport

## Which information?

Method to isolate the influence of mass transfer from the electrochemical impedance response of a system

## What was the aim of this overview?

EIS is only a part of what is possible to do with Function Transfer.  
Other methods exist and maybe others can be developed to investigate the dynamic behavior of the system.

User has to find the appropriate input and output.

	System	Information
<b>ElectroHydrodynamic TF (EHD)</b>	<ul style="list-style-type: none"> <li>• Electrochemical system that are either partially or completely limited by mass transport</li> <li>• Partially blocked electrode</li> </ul>	<ul style="list-style-type: none"> <li>• Mass transfer from the electrochemical impedance response of a system</li> </ul>
<b>PhotoElectrochemical TF (IMVS/IMPS)</b>	<ul style="list-style-type: none"> <li>• Solar cell</li> </ul>	<ul style="list-style-type: none"> <li>• Electron lifetime</li> <li>• Diffusion</li> </ul>
<b>MagnetoHydrodynamic TF</b>	<ul style="list-style-type: none"> <li>• Metal electrode deposition</li> </ul>	<ul style="list-style-type: none"> <li>• Kinetic of the processes</li> </ul>
<b>ElectroGravimetric TF</b>	<ul style="list-style-type: none"> <li>• Insertion in film</li> <li>• Reaction with adsorbed species</li> </ul>	<ul style="list-style-type: none"> <li>• Kinetic of the processes</li> <li>• Chemical identification of the species</li> </ul>
<b>Raman TF</b>	<ul style="list-style-type: none"> <li>• Reaction with adsorbed species</li> </ul>	<ul style="list-style-type: none"> <li>• dynamic information on the interface</li> </ul>
<b>ThermoElectrochemical TF (TEC)</b>	<ul style="list-style-type: none"> <li>• Redox system</li> </ul>	<ul style="list-style-type: none"> <li>• kinetic of mass transport</li> </ul>
<b>PneumatoChemical Impedance Spectroscopy</b>	<ul style="list-style-type: none"> <li>• Insertion reaction</li> <li>• Hydrogen insertion</li> </ul>	<ul style="list-style-type: none"> <li>• phase transformation</li> </ul>

Feel free to visit our web site, some application notes or EIS handbook may be helpful for your applications:

<http://www.bio-logic.info/potentiostat/notesan.html>

Thank you for your attention

Lets move to the instruments