

80 A dedicated to electrochemical applications requiring current

- FUEL CELLS
- SUPERCAPACITORS
- ELECTROPLATING
- ELECTROLYSIS

With the increasing investment in new energy sources, researchers need high current test stations. The **HCP-803** addresses these applications requiring high currents and also general electrochemistry applications. **HCP-803** is designed for fuel cell aging tests (single cell or small stacks), for the characterization and performances study of supercapacitors and for electrochemical plating experiments.



This instrument is a research-grade potentiostat with a 80 A current booster option. It is available in two different models: the controlled stand-alone **HCP-803** and an external booster **VMP3B-80** that is included in another chassis to be controlled by an existing potentiostat of our range. This booster is compatible with **SP-150**, **VSP** and **VMP3**.

Our external current booster is plug-and-play. It can be connected and disconnected from the channel board and reconnected to another one without switching off the instrument.

You can connect the **HCP-803** directly to the PC with a USB port 2.0. The alternative Ethernet communication allows installation on a Local Area Network. Analog inputs/outputs are available to manage external instruments and record the generated data.

The **EC-Lab**[®] software supplied with the unit offers more than 50 techniques that can be sequenced or linked, and also a variety of analysis tools.

EIS capability is provided in standard with dedicated techniques to control the cell in potentiostatic or galvanostatic mode. Among the different software analyses, modelling tool with Levenberg-Marquardt and Simplex algorithms is interesting to study materials constituting the cells and for ageing follow up.

GENERAL SPECIFICATIONS

- High current system: ± 80 A
- Voltage range: ± 3 V @ ± 80 A
and ± 4 V @ ± 40 A
- EIS capability from 10 μ Hz to 10 kHz and up to 1 MHz without booster
- EC-Lab[®] and EC-Lab[®] Express software
- Disconnectable current booster to use the system as a research grade potentiostat/galvanostat
- Available as a 80 A/3 V booster, compatible with SP-150, VSP and VMP3

EC-LAB® EXPRESS: EASY-TO-USE AND POWERFUL MONITORING SOFTWARE

EC-Lab® Express software offers a full range of electrochemical techniques (more than 30 linkable techniques containing up to 100 sequences).


This software is very easy to use. The settings and graphs are shown on one screen view. An experiment selector enables the user to quickly switch between techniques. The user can set data sampling and recording conditions without any limit on the number of data points taken.



EC-LAB® A COMPLETE SOFTWARE FOR FULL CONTROL OF THE EXPERIMENT

Especially devoted to energy devices testing, EC-Lab® software offers more than 50 techniques with up to 100 sequences that can be linked.

Many experimental parameters can be modified "on the fly" during the experiment, with the change stored into the raw data file. The appearance of the software interface is able to be adapted to create the best working environment for the user. A powerful technique builder can execute a series of different modular techniques and loop options to create complex experimental sequences.



In case of fuel cell testing, the monitored current can be used to control gas flow meters thanks to the external device configuration window.

ANALYSIS TOOLS

The complete graphic package provided with EC-Lab® software includes analysis tools and advanced fitting tools.

A "process data" function helps the user to calculate additional variables such as energy, charge or capacity during successive cycles.

EC-Lab® software offers classical analysis tools (linear or circular fit) and also a powerful tool for EIS data fitting. It includes electrical equivalent circuits with basic electronic elements and uses two minimization algorithms (Simplex, Levenberg-Marquardt).

TECHNIQUES

Voltammetric techniques:
OCV, CV, CVA, CA, CP, SV

Impedance Spectroscopy:
Galvano/potenti EIS, staircase EIS (Mott-Schottky), multisine mode

Technique builder:
Modular potenti/galvano, loop, trigger in/out, wait

Batteries and Supercapacitors:
Galvano/potenti cycling, constant load/power discharge, profile simulation

Fuel cell testing:
I-V Characterisation, constant load, constant power

Ohmic Drop:
Manual IR, EIS determination

Stack techniques:
I-V characterization, potenti/galvano EIS, constant current, constant voltage

SPECIFICATIONS

CHANNEL BOARD

Cell control

Connection	2, 3, 4 or 5 terminal leads (+ ground)
Compliance	10 V range adjustable from ± 10 V to 0 - 20 V
Maximum current	± 400 mA continuous
Maximum potential resolution	300 μ V on 20 V dynamic range programmable down to 5 μ V on 200 mV range
Maximum current resolution	0.004 % of the dynamic range programmable down to 760 pA on the 10 μ A range
Accuracy (DC)	< 0.1 % FSR*
Rise time	< 2 μ s (no load)
Acquisition time	20 μ s

Current measurement

Ranges	± 10 μ A to ± 400 mA (7 ranges)
Maximum resolution	0.004 % FSR*
Acquisition speed	200,000 samples/s
Accuracy (DC)	< 0.1 % FSR*

Potential measurement

Ranges	± 2.5 V, ± 5 V, ± 10 V, ± 10 V adjustable
Maximum resolution	0.0015 % of the range, down to 75 μ V
Acquisition speed	200,000 samples/s
Accuracy (DC)	< 0.1 % FSR*

Electrometer

Inputs	3 potential measurements
Impedance	> 10^{12} ohms in parallel with < 20 pF
Bias current	< 5 pA

Auxiliary inputs/outputs

2 analog inputs	automatic ± 2.5 V, ± 5 V, ± 10 V ranges
1 analog output	± 10 V
1 input trigger	TTL level
1 output trigger	TTL level

General

Dimensions, weight	260 x 495 x 465 (mm, H x W x D), 23 kg
Power	200-240 V, 47-440 Hz

IMPEDANCE (EIS)

Impedance

Frequency range	10 μ Hz to 10 kHz (see contour map)
Amplitude	1 mVpp to 1 Vpp
	0.1 % to 50 % of the current range
Accuracy	2 %, 2°

BOOSTER BOARD

Cell control

Connection	5 terminal leads
Compliance	-5; +5 V range (see contour map)
Maximum current	80 A continuous (see contour map)
Potential ranges	± 5 V @ 1 A, ± 4 V @ 40 A and ± 3 V @ 80 A
Rise time Potentio	95 μ s
Galvano	150 μ s

Measurement

Potential accuracy (DC)	< 0.1% FSR*
Current accuracy (DC)	< 0.5 % FSR*
Current noise (peak to peak 0-100 kHz)	14 mArms
Potential noise (peak to peak 0-100 kHz)	0.18 mVrms


Electrometer

Inputs	3 potential measurements
Impedance	10^{10} Ohms

Auxiliary inputs/outputs

1 security input to open circuit
1 emergency stop push button

CONTOUR MAP



* FSR: Full Scale Range

Specifications subject to change

Bio-Logic, SAS
1, rue de l'Europe
38640 CLAIX - France
Tel.: +33 476 98 68 31
Fax: +33 476 98 69 09
www.bio-logic.info



Bio-Logic USA, LLC
P.O.Box 30009
Knoxville TN 37930 - USA
Tel: +1 865 769 3800
Fax: +1 865 769 3801
www.bio-logic.us