

Dynamic resistance determination: Is there a link between EIS and DC measurements?

Aymeric Pellissier, Jean-Paul Diard, Sébastien Benoit,
Bio-Logic Science Instruments

INTRODUCTION

Currently, several types of batteries are investigated (Li-ion, LiFePO₄, ...). Each type has its own specifications such as nominal voltage (V), capacity (A.h), power (W), energy (W.h), cyclability... The dynamic resistance R_{dt} of a battery is also a relevant parameter. This resistance can be determined by a simple current interrupt technique. In this technique, a current pulse is applied to the battery and the voltage shift is measured after a certain period of time Δt (Fig. 1).

The ratio of the voltage shift to the current difference ΔI yields to a resistance, R_{dt} :

$$R_{dt} = \Delta E(t) / \Delta I(t) \quad [1]$$

The duration Δt may be different depending on the used procedure and the instrument. For example, this period is 30 s in the procedure of the United States Advanced Battery Consortium (USABC).

The aim of this poster is to compare the resistance measured by Electrochemical Impedance Spectroscopy (EIS) and DC measurements.




Fig. 1: Current and voltage change versus time after current pulse applied at t₁. ΔE(t) = E₂-E₁; ΔI(t) = I₂-I₁

DYNAMIC RESISTANCE DETERMINATION...

During a discharge (Fig. 2), the dynamic resistance is determined by DC method and EIS measurement. The DC and EIS measurements are shown in Figs. 3 and 4, respectively




Fig. 2: Voltage evolution versus charge. The EIS and DC measurements are indicated by the red arrows. (Inset: Voltage vs. time).




Fig. 3: Dynamic resistance change at certain period after the current pulse: Δt: 200 μs, 1 ms; 1.6 s and 10 s.




Fig. 4: Nyquist plot (inset: Bode plot).

... AND COMPARISON

To compare the data obtained with DC or EIS experiments, the dynamic resistance and impedance after each discharge step are plotted versus time (Fig. 5). Both DC and AC data have to be given in terms of frequency.

The DC data conversion is done by $\omega = 1/\Delta t$ where ω is the pulsation in rad.s⁻¹.

So if $\Delta t = 200 \mu s$, the corresponding frequency is 5 kHz.

It is noteworthy that the impedance is not a pure resistance at 5 kHz (not only real part). The phase at this frequency is -20° (Fig. 4).

- similar trends (small increase)
- same magnitude (between 25-50 mΩm).

FURTHER INTO THE THEORY

In order to compare the time domain and the frequency domain measurement, it is possible to calculate the step response of the battery using the inverse of the Laplace Transform, L⁻¹.

For example, let us consider the equivalent circuit (typical for a battery):

$$Z(s) = R_1 + L_1 s + R_2 / (1 + R_2 C_2 s) + \sigma / s^{1/2} \quad [2]$$

For a step of current, ΔI, the step potential response of the circuit is given by:

$$E(t) = L^{-1}[(\Delta I/s)Z(s)] = \Delta I (R_1 + L_1 \delta(t) + R_2 [1 - \exp(-t/(R_2 C_2)) + 2\sigma t^{1/2}/\pi^{1/2}]) \quad [3]$$

Where δ(t) is the Dirac delta function.

The ratio E(t)/ΔI(t) and modulus are given by the Eqs. (3) and (4), respectively.

$$|Z(\omega)| = [(R_1 + \sigma/\omega^{1/2} + R_2/(1+R_2^2 C_2^2 \omega^2))^2 + (R_1 + \sigma/\omega^{1/2} + R_2/(1+R_2^2 C_2^2 \omega^2))^2]^{1/2} \quad [4]$$

Modulus versus angular frequency (from AC measurement) and ratio E(t)/ΔI(t) versus the inverse of Δt (from DC measurement) are plotted in Fig. 6. The superimposition of both plots show that data coming from AC and DC measurement are in good agreement.

CONCLUSION

The dynamic resistance determined by DC method, for instance GITT techniques, can be compared to impedance data. The conclusion of this comparison is that the dynamic resistance can be used as a rough approximation of the impedance measured by AC method. The values of the dynamic resistance or the modulus are theoretically linked together via Eq. (3) and lead to similar trends.

Nonetheless, to fully characterize the battery, an investigation on the whole frequency range by AC method is needed.

More info on: <http://www.bio-logic.info/potentiostat/notesan.html>




Fig. 5: Dynamic resistance and impedance plots resulting from DC (200 μs) and Dynamic resistance AC (~5 kHz) investigations, respectively.




Fig. 6: |Z| and Dynamic resistance (DC simulation) plot versus logarithm of ω or Δt, respectively (inset: The corresponding Nyquist plot).