

# CASP, VASP : NEW NON-LINEAR IMPEDANCE TECHNIQUES FOR THE DETERMINATION OF GENERALIZED CORROSION PARAMETERS

N. PORTAIL <sup>(1)</sup>, S. BENOIT <sup>(1)</sup>, J.-P. DIARD <sup>(1,2)</sup>, C. MONTELLA <sup>(2)</sup>

1. Bio-Logic SAS, Claix

2. Laboratoire d'Electrochimie et de Physico-chimie des Matériaux et des Interfaces, Grenoble

## INTRODUCTION

When a metal is corroding under the action of an oxidizing species without mass transport limitation (*i.e.* Tafelian system) the steady-state  $I$  vs.  $E$  curve follows the Stern or Wagner-Traud relationship :

$$I = I_{corr} \left( \exp\left(\frac{E - E_{corr}}{\beta_a / \ln(10)}\right) - \exp\left(\frac{-(E - E_{corr})}{\beta_c / \ln(10)}\right) \right) \quad (1)$$

where  $E_{corr}$  and  $I_{corr}$  are the corrosion potential and current, respectively,  $\beta_a$ ,  $\beta_c$ , the Tafel parameters. The classical methods (Tafel,  $R_p$ ) based on the use of the equation (1) give an estimation of the generalized corrosion current of a metal immersed in an aggressive solution. Two new methods are available in EC-Lab<sup>®</sup>, that are founded on non-linear impedance measurements. The first one (CASP) analyzes the current response harmonics and the second one (VASP) the evolution of the polarization resistance  $R_p$  with the amplitude of a sinusoidal potential modulation around  $E_{corr}$ .

## METHODS / DESCRIPTION

### CASP: Constant Amplitude Sinusoidal microPolarization

$$|i(t)| = I_{corr} \left( \exp\left(\frac{V_a \sin(2\pi f_s t)}{\beta_a / \ln(10)}\right) - \exp\left(\frac{-V_a \sin(2\pi f_s t)}{\beta_c / \ln(10)}\right) \right) \quad (2)$$

- The current response is a periodic signal, expandable as a Fourier series.
- The number of apparent harmonics depends on the amplitude  $V_a$  of the potential modulation.
- Only the three first harmonics  $\delta I_1$ ,  $\delta I_2$ ,  $\delta I_3$  are used.

$$I_{corr} = \frac{1}{4\sqrt{3}} \frac{|\delta I_1 + 3\delta I_3|^2}{\sqrt{|\delta I_2^2 + 2\delta I_3| |\delta I_1 + 3\delta I_3|}} \quad (3)$$

$$b_{a,c} = \frac{2}{V_a} \frac{\sqrt{3} \sqrt{|\delta I_2^2 + 2\delta I_3| |\delta I_1 + 3\delta I_3|} \pm \delta I_2}{\delta I_1 + 3\delta I_3} \quad (4)$$

## RESULTS / DISCUSSION

Ni in 0,5 mol L<sup>-1</sup> H<sub>2</sub>SO<sub>4</sub>




Fig. 3 : Top : Electrode potential modulation, Bottom : Current response (Ni in 0.5 mol L<sup>-1</sup> H<sub>2</sub>SO<sub>4</sub>)




Fig. 4 : Discrete Fourier Transform of the current response Fig. 3  $V_a = 25$  mV

The impedance of the system at the frequency measurement must be close to  $R_p$ .

Using relationships (3, 4), it is possible to obtain the Tafel parameters:  $I_{corr} = 10,6 \mu A$ ,  $\beta_a = 84,05$  mV,  $\beta_c = 353,3$  mV.

## CONCLUSIONS

The CASP and VASP methods can be used to measure the corrosion current of a metallic electrode undergoing uniform corrosion. These methods carry out impedance measurements for which the amplitudes of the modulation correspond to a non-linear behaviour. Thanks to these methods, the determination of the Tafel parameters can be avoided in the calculation of  $I_{corr}$  (using the Stern-Geary equation).

These methods can be extensively used for the monitoring of corroding materials.

More info can be found in the corresponding application notes:  
<http://www.bio-logic.info/potentiostat/notesan.html>

### VASP: Variable Amplitude Sinusoidal microPolarization

- It is based on the measurement of the evolution of  $R_p$  with the amplitude of the potential modulation  $V_a$ .
- $R_p$  is determined by EIS (Electrochemical Impedance Spectroscopy) at a low enough fixed frequency  $f_s$ .

$$\frac{1}{R_n} = I_{corr} \sum_{k=0}^{\infty} \frac{b_a^{2k+1} + b_c^{2k+1}}{2^{2k} k! |k+1|!} V_a^{2k} \quad (5)$$




Fig. 1 : Sine train with an amplitude varying from  $V_{amin}$  to  $V_{amax}$  with a constant frequency  $f_s$




Fig. 2 : Evolution of the impedance of a non-linear system (obtained on test box 3#2) with an increasing amplitude  $V_a$

Ni in 0,1 mol L<sup>-1</sup> HCl




Fig. 5 :  $R_p$  vs.  $V_a$  (Ni in 0,1 mol L<sup>-1</sup> HCl)

For a negligible electrolyte resistance value  $R_{\Omega}$ , the polarization resistance is given by equation (5).

Tafel parameters are determined by a parametrical identification:  $I_{corr} = 11,7 \mu A$ ,  $\beta_a = 74,7$  mV,  $\beta_c = 190,4$  mV.